

# Nexfin® non-invasive continuous hemodynamic monitoring: validation against continuous pulse contour and intermittent transpulmonary thermodilution derived cardiac output in critically ill patients

Koen Ameloot, Katrijn Van De Vijver, Ole Broch, Niels Van Regenmortel, Inneke De Laet, Karen Schoonheydt, Hilde Dits, Berthold Bein, Manu L. N. G. Malbrain

K. Ameloot, K. Van De Vijver, N. Van Regenmortel, I. De Laet, K. Schoonheydt, H. Dits, M. L. N. G. Malbrain (✉) Department of Intensive Care Medicine. Ziekenhuis Netwerk Antwerpen. ZNA Stuivenberg, Lange Beeldekensstraat 267, B-2060 Antwerpen 6, Belgium, Tel: +32 3 / 217 7399, Fax: +32 3 /217 7574, manu.malbrain@skynet.be

B. Bein, O. Broch, Department of Anaesthesiology and Intensive Care Medicine, University Hospital Schleswig-Holstein, Campus Kiel, Schwannenweg 21, D-24105 Kiel, Germany, Tel.: +49 431 / 597 2971, Fax: +49 431 / 597 3002, bein@anaesthesie.uni-kiel.de

Received: 14.09.2013, accepted: 22.10.2013

**Abstract Introduction** Nexfin® (Bmeye, Amsterdam, The Netherlands) is a completely non-invasive blood pressure and cardiac output (CO) monitor based on finger arterial pulse contour analysis. The aim of this study was to validate this technique in a mixed population of critically ill patients. **Patients and methods** We performed a prospective open observational study in a mix of medical/surgical and burns critically ill patients (n=45) to validate Nexfin® against transpulmonary thermodilution and pulse contour CO (PiCCO, Pulsion Medical Systems, Munich, Germany). Nexfin® cardiac output (NexCO) and PiCCO pulse contour CO (CCO) were measured continuously and recorded at 2 hour intervals during the 8-hours study period. PiCCO thermodilution CO (TDCO) was measured at 0—4—8 hours. Statistical analysis was performed by Pearson regression, Bland and Altman, Concordance plot and Polar plot analysis. **Results** NexCO showed a moderate to good (significant) correlation with TDCO ( $R^2=0.68, p<0.001$ ) and CCO ( $R^2=0.71, p<0.001$ ). Bland and Altman analysis comparing NexCO with TDCO revealed a bias ( $\pm$ limits of agreement, LA) of  $0.4\pm 2.32$  L/min (with 36% error) while analysis comparing NexCO with CCO showed a bias ( $\pm$ LA) of  $0.2\pm 2.32$  L/min (37% error). Subgroup analysis showed that NexCO keeps reasonable performance even in unstable septic patients with a high CO, low SVRI or on high dose of norepinephrine. NexCO is able to follow relative changes in TDCO and CCO during the same time interval (level of concordance 89.3% and 81% respectively). The absolute amplitude correlation of these changes was clinically insufficient but statistically significant ( $R^2=0.63, p<0.001$  and  $R^2=0.31, p=0.006$  respectively). Finally, polar plot analysis showed that trending capabilities were acceptable when changes in NexCO ( $\Delta$ NexCO) were compared to  $\Delta$ TDCO and  $\Delta$ CCO (respectively 89% and 88.9% of changes were within the level of 10% limits of agreement). Conclusion in this sample of critically ill patients we found a moderate to good correlation between CO measurements obtained with Nexfin® and PiCCO. The Nexfin® can be used for keeping track of changes in CO over time (for instance to assess the therapeutic effect of a given treatment), although the absolute criteria for full interchangeability were not met.

**Key words** cardiac output · hemodynamic monitoring · minimal invasive · Nexfin® · PiCCO · pulse contour analysis · thermodilution

## Introduction

The true value of continuous hemodynamic monitoring in critically ill patients becomes clear in the light of beat-to-beat changing hemodynamics due to either continuously improving or deteriorating cardiac and disease status [6, 16]. Therefore, the

critically ill patient must be resuscitated to a continuously changing optimal left ventricular end-diastolic volume. This should be titrated together with an accurate dose of vasopressor agents and inotropes to optimize circulation and restore end-organ perfusion without causing harm due to excessive fluids [20]. Also, for further decision-making we often want

to measure in real-time the hemodynamic effects of ongoing therapeutic interventions. The PiCCO system (PiCCO2, Pulsion Medical Systems, Munich, Germany) uses a dedicated PiCCO thermistor-tipped arterial catheter to analyze the patient's heart rate and arterial pressure waveform continuously [24]. Due to the unique properties of each patient's arterial tree, initial calibration of the monitoring system using transpulmonary thermodilution CO measurement (TDCO) improves accuracy of the beat-to-beat cardiac output (CO) obtained by pulse contour analysis (CCO) [13]. The PiCCO device has been validated in numerous studies including burns, medical and surgical critically ill patients [3, 15–17, 21, 28, 30]. However, PiCCO remains a relative invasive technique that requires both an arterial and a central venous catheter and therefore increases the risk of iatrogenic complications such as pneumothorax, bleeding, catheter sepsis and deep venous thrombosis [1]. Moreover, the need for calibration by transpulmonary thermodilution may delay the initial measurement and is time consuming, the system is cost intensive and cannot be used pre-hospital or in a regular ward. The Nexfin® (BMEYE, Amsterdam, the Netherlands) device is a totally non-invasive continuous blood pressure and CO (NexCO) monitor based on finger arterial pressure pulse contour analysis. Potential advantages include its non-invasiveness and ease of use. Questions rose on the accuracy of Nexfin® measurements in critically ill, as there is no initial calibration of the monitoring system to adjust for the unique mechanical properties of each patient's arterial tree. Moreover, there are concerns about the reliability of continuous non-invasive finger blood pressure derived pulse contour analysis in patients with reduced perfusion of the hand due to high systemic vascular resistance (SVR) or hypothermia. We performed an open observational study in a mix of medical/surgical and burns critically ill patients to validate Nexfin® against both TDCO and CCO obtained by PiCCO.

## Methods

### Study Population

We prospectively studied 47 critically ill patients admitted to the intensive care units of our hospital. Inclusion criteria were hemodynamic instability with need for continuous hemodynamic monitoring and the presence of a central venous (jugular or subclavian) catheter and a PiCCO femoral arterial catheter already in place before inclusion in the study.

### Nexfin® technique

The Nexfin® (BMEYE, Amsterdam, The Netherlands) method is based on the measurement of finger arterial pressure by an inflatable cuff around the middle phalange of the finger. The pulsating finger artery is clamped to a constant volume by applying a varying

counter pressure equivalent to the arterial pressure using a built-in photoelectric plethysmograph and an automatic algorithm (Physiocal). The resulting finger arterial pressure waveform is reconstructed into a brachial artery pressure waveform by a generalized algorithm. NexCO is calculated by a pulse contour method (CO-TREK) using the measured systolic pressure time integral and the heart's afterload determined from the Windkessel model [33].

## Measurements

In an 8-hours period, simultaneous CCO and NexCO measurements were obtained every 2 hours (0–2–4–6–8 hours in total 225 paired measurements) while simultaneous TDCO and NexCO were obtained every 4 hours (0–4–8 hours, in total 135 paired measurements). The CCO and NexCO values were recorded simultaneously by hand 5 minutes before TDCO was determined by 3 repeated injections of 20 ml of sterile ice-cold saline via the central venous line. Blood pressure measurements were recorded continuously by Nexfin® and PiCCO and were by each device used to calculate the continuous CO. Sub-analysis was performed for patients with a low MAP (defined as  $\text{MAP} \leq 70 \text{ mmHg}$ ), low and high TDCO (defined as  $\leq 4 \text{ L/min}$  and  $\geq 8 \text{ L/min}$ ), low and high SVRI (defined as an SVRI obtained by PiCCO  $\leq 1700 \text{ dyne.s.cm}^{-5}/\text{m}^2$  and  $\geq 3000 \text{ dyne.s.cm}^{-5}/\text{m}^2$ ), patients on high dose norepinephrine ( $\geq 0.3 \mu\text{g/kg/min}$ ) and hypothermia ( $T \leq 35^\circ\text{C}$ ). To assess the ease of use we measured the time to initial measurement, the number of repositions needed in the 8 hours observation period and also the nurses filled in a questionnaire (see addendum) ( $n=27$  patients).

## Data analysis and statistics

Results are presented as mean ( $\pm$ SD) unless otherwise stated. Mean values were compared using student's t-test. Paired CO measurements by 2 different methods were compared statistically using 4 different methods. First, we used Pearson correlation and linear regression analysis. Two methods are considered equal if the line of identity crosses the origin of X and Y-axis and if  $R^2$  ( $R$ =Pearson's correlation coefficient) is  $>0.6$ . Second, we calculated bias, precision and limits of agreement (Bland-Altman analysis [4]) and the percentage error (PE, defined as two times SD of the bias over the mean TDCO or CCO) as described by Critchley and Critchley [10]. If the differences within bias  $\pm 1.96 \text{ SD}$  (limits of agreement, LA) are not clinically important, if the precision of the new technique is comparable to the reference technique and if the percentage error is less than 30% the two methods may be used interchangeably [7]. Third, the ability of NexCO to track changes or trends in TDCO or CCO was assessed by plotting  $\Delta\text{TDCO}$  or  $\Delta\text{CCO}$  against  $\Delta\text{NexCO}$  during the same time interval (four quadrants trend plot). The concordance is calculated as the percentage of pairs with the same direction of




Fig. 1. Cardiac output measurements: TDCO vs NexCO. Only one average value per patient is plotted. Panel A. Regression analysis. Panel B. Bland-Altman analysis. Patient averages with the mean cardiac output ranges (X-axis) and errors (Y-axis) during the 8 hour study period. Dotted line indicates mean bias and solid lines indicate lower and upper limit of agreement. NexCO: Nexfin® cardiac output, TDCO: thermodilution cardiac output.




Fig. 2. Cardiac output measurements: CCO vs NexCO. Only one average value per patient is plotted. Panel A. Regression analysis. Panel B. Bland-Altman analysis. Patient averages with the mean cardiac output ranges (X-axis) and errors (Y-axis) during the 8 hours study period. Dotted line indicates mean bias and solid lines indicate lower and upper limit of agreement. CCO: pulse contour continuous cardiac output, NexCO: Nexfin® cardiac output.

change. Based on previous reports, the concordance should be >90% when pairs with both a  $\Delta$ TDCO or  $\Delta$ CCO and  $\Delta$ NexCO  $\leq \pm 1$  L/min (or less than 15% of change) are excluded for analysis [2]. Finally trending capability of the NexCO compared to TDCO and CCO was assessed by polar plots as suggested by Critchley et al. [11]. Concordance analysis looking at direction of changes is a very simple but crude measure of how well 2 measurements trend. Important aspects of the measurement, such as the magnitude of the underlying CO change ( $\Delta$ CO), and the degree of agreement, are totally ignored. Therefore, Critchley et al. suggested converting the x-y values to polar coordinates, where agreement is shown by the angle the vector makes with the line of identity ( $y=x$ ) and magnitude of change by the length of the vector [11]. Thus, statistical measures that fully represent the magnitude of  $\Delta$ CO and its degree of agreement are retained. From these data, a new polar plot that shows agreement as the angle  $\Delta$  (angle made by  $\Delta$ CO vector with the line of identity [ $y=x$ ]) against the change in CO as the radian (distance of data point from center of polar plot) can be drawn. Conversion of the Cartesian data with regard to change in cardiac output ( $\Delta$ CO) into a

(x, y) polar coordinate format was performed using an Excel spreadsheet (Microsoft Office Excel 2007; Microsoft Corp.), with the formulas for calculation of absolute value of mean  $\Delta$ CO, quadrant, radians and angle as suggested by Critchley et al. [11].

### Ethics

The study was conducted in accordance with the ICU protocol, the declaration of Helsinki and applicable regulatory requirements as approved by the institutional review board and the local institutional ethics committee (approval number 3789). In view of the nature of the study being purely observational, not demanding a deviation from standard clinical ICU care and since the results obtained by Nexfin® were not used for clinical decision making, informed consent from the patient or the next of kin was not deemed essential. We merely analysed the existing situation and did nothing to influence events. Only treating ICU physicians accessed the medical records. All data were pseudonymized before analysis.




Fig. 3. Four quadrants trend plot. Panel A. Plot for 90 paired measurements of  $\Delta$ NexCO and  $\Delta$ TDCO. From the 90 initial paired measurements, 34 pairs were excluded (exclusion zone is indicated as grey dots within grey-shaded square) because either  $\Delta$ NexCO and  $\Delta$ TDCO were  $\leq \pm 15\%$  or because  $\Delta$ NexCO or  $\Delta$ TDCO was equal to zero. The calculated level of concordance was 89.3% (50/56) (6 pairs fell within the upper left or lower right quadrant and correspond to poor concordance, black dots). See text for explanation. Panel B. Plot for 180 paired measurements of  $\Delta$ NexCO and  $\Delta$ CCO. From the 180 initial paired measurements, 75 pairs were excluded (exclusion zone is indicated as grey-shaded square) because either  $\Delta$ NexCO and  $\Delta$ CCO were  $\leq \pm 15\%$  change or because  $\Delta$ NexCO or  $\Delta$ CCO was equal to zero. The calculated level of concordance was 81% (85/105). See text for explanation.

## Results

### Study population

In 2 patients (4.3%) it was impossible to obtain Nexfin<sup>®</sup> values from any of 10 fingers and they were therefore excluded from final analysis. Baseline characteristics of the 45 remaining patients are summarized in table 1. Thirty-one patients (69%) were in shock (reflected by an elevated arterial lactate) with the majority in septic shock (n=18, 58%). Only a minority was in cardiogenic shock (n=6, 19%). A total of 35 patients (78%) received norepinephrine at a mean ( $\pm$ SD) dose of  $0.20 \pm 0.17$  (range: 0.02—1)  $\mu\text{g}/\text{kg}/\text{min}$ , while 27 patients (60%) received dobutamine at a dose of  $4.30 \pm 2.10$  (range: 1—10)  $\mu\text{g}/\text{kg}/\text{min}$  (range: 1—10). Forty-three patients (96%) were mechanically ventilated and the 2 remaining patients were non-invasively ventilated. The critically illness of the patient sample is reflected by high scores on 3 different ICU scoring systems (APACHE II, SOFA and SAPS II). In-hospital mortality was 57.8%. The neurological, respiratory and hemodynamic parameters and the dose of the infused drugs did not show significant changes during the entire study period (Table 2).

### Cardiac Output

Mean NexCO was comparable to mean TDCO ( $6.1 \pm 2.3$  vs  $6.6 \pm 2.2$  L/min,  $p=0.10$ ) and to mean CCO ( $6.1 \pm 2.3$  vs  $6.4 \pm 2.3$  L/min,  $p=0.30$ ). Correlation, regression and the Bland and Altman analysis are shown in Table 3, Table 4 and Figure 1 and Figure 2. Pearson correlation coefficients comparing NexCO with TDCO ( $R^2=0.68$ ) and NexCO with CCO ( $R^2=0.71$ ) were comparable and showed a highly significant (both  $p<0.001$ ) correlation between all obtained CO measurements. Bland and Altman analysis com-


paring NexCO with TDCO revealed a mean bias $\pm$ LA of  $0.4 \pm 2.32$  L/min (with 36% error) while analysis comparing NexCO with CCO showed a bias ( $\pm$ LA) of  $0.2 \pm 2.32$  L/min (37% error). TDCO was highly correlated with CCO ( $R^2=0.95$ ,  $p<0.001$ ) with a bias ( $\pm$ LA) of  $0.2 \pm 0.86$  L/min (13.3% error).

Subanalysis for patients with a low MAP, high TDCO, low SVRI and high dose norepinephrine consistently showed a very good correlation between NexCO and TDCO or CCO. NexCO was less reliable in patients with hypothermia and not reliable in patients with low TDCO and high SVRI (Table 3, Table 4).

The four quadrants concordance plots are shown in Figure 3. From the 90 initial paired  $\Delta$ NexCO/ $\Delta$ TDCO measurements, 34 pairs were excluded because either  $\Delta$ NexCO or  $\Delta$ TDCO were  $\leq \pm 1\text{L}/\text{min}$  (or  $\leq 15\%$  change) or because  $\Delta$ NexCO or  $\Delta$ TDCO were equal to zero (panel A). The calculated level of concordance was 89.3% (50/56). The absolute amplitude correlation of these changes was clinically sufficient ( $R^2=0.63$ ,  $p<0.001$ ). From the 180 initial paired  $\Delta$ NexCO/ $\Delta$ CCO measurements, 75 pairs were excluded because either  $\Delta$ NexCO or  $\Delta$ CCO were  $\leq \pm 1\text{L}/\text{min}$  (or  $\leq 15\%$  change) or because  $\Delta$ NexCO or  $\Delta$ CCO were equal to zero (panel B). The calculated level of concordance was only 81% (85/105). The absolute amplitude correlation of these changes was clinically insufficient but still significant ( $R^2=0.31$ ,  $p=0.006$ ).

The polar trending plots are shown in Figure 4. From the 90 initial data 98.9% of the data points were within the 20% lines and 89% within the 10% lines, suggesting acceptable trending capabilities (panel A). From the 180 initial data 98.3% of the data points were within the 20% lines and 88.9% within the 10% lines, suggesting acceptable trending capabilities (panel B).

Panel A



Panel B




Fig. 4. Polar plot. The distance from the center of the plot represents the mean change in cardiac output ( $\Delta\text{CO}$ , expressed as %, with 1,0 referring to 100% change from baseline) and the angle  $\Delta$  with the horizontal (0-degree radial) axis represents agreement. The less the disagreement between CO measurements, the closer data pairs will lie along the horizontal radial axis. Data with good trending will lie within 10% limits of agreement. However, data with poor trending will be scattered throughout the plot and lie outside the limits of good and acceptable agreement (i.e., 10% and 20%, respectively). See text for explanation. Panel A. Polar plot for 90 paired measurements of mean  $\Delta\text{CO}$  (%), calculated as absolute value of  $(\Delta\text{NeXCO} + \Delta\text{TDCO})/2$ . From the 90 initial data 98.9% of the data points lie within the 20% lines and 89% within the 10% lines, suggesting acceptable trending capabilities. Panel B. Polar plot for 180 paired measurements of mean  $\Delta\text{CO}$  (%), calculated as absolute value of  $(\Delta\text{NeXCO} + \Delta\text{CCO})/2$ . From the 180 initial data 98.3% of the data points lie within the 20% lines and 88.9% within the 10% lines, suggesting acceptable trending capabilities.

Table 1. Patient characteristics.

	Mean	$\pm$	SD	n=45 (100%)
Demographics				
Age (yrs)	57.6	$\pm$	19.4	
Male				32 (71%)
Reason of admission				
Medical				27 (60%)
Surgical				9 (20%)
Trauma				5 (11%)
Burns				4 (9%)
Shock				31 (69%)
Septic				18 (58%)
Cardiogenic				6 (19%)
Other				7 (23%)
ICU scores				
APACHE II	25.3	$\pm$	10.3	
SOFA	9.4	$\pm$	3.3	
SAPS II	51.5	$\pm$	16.9	

#### Ease of use

Data on ease of use were collected in 27 patients. There were no local signs of disturbed circulation in the middle finger due to the application of the finger cuff. The time between the decision to apply Nexfin<sup>®</sup> and the first measurement was less than 5 minutes in 23/27 patients (85%) and between 5–10 minutes

in 4/27 patients (15%). In 9/27 patients (33%) we were able to do measurements with the first application, while 13/27 patients (48%) needed 1–5 and 5/27 (19%) needed more than 5 repositions. Nurse questionnaires revealed a mean score of  $1.4 \pm 0.5$  for the set-up of the device,  $1.7 \pm 0.7$  for set-up placement,  $1.8 \pm 0.5$  for measurements and  $1.9 \pm 0.5$  for ease of use (1=very easy to 4=very difficult).

#### Discussion

We performed an open observational study in 45 mixed surgical/medical and burns critically ill patients to validate the Nexfin<sup>®</sup> against transpulmonary thermodilution and continuous femoral arterial pulse contour derived CO by the PiCCO. To the best of our knowledge, this is the second Nexfin<sup>®</sup> CO validation study conducted in mainly medical ICU patients, the first being published last year by Monnet and co-workers [22].

First, we found moderate to good CO correlation coefficients with TDCO ( $R^2=0.68$ ) and CCO ( $R^2=0.71$ ). The overall calculated PE's were however too high to meet the criteria for general interchangeability as suggested by Critchley et al [11]. These results are in line with previous Nexfin<sup>®</sup> CO validation studies against TEE and PiCCO during abdominal and cardiac surgery [5, 9] and against PAC in a small sample of 10 post-surgical ICU patients [31]. In 2 studies a  $PE < 30\%$  was reported [5, 31]. Our PE in critically ill patients was lower than the one reported previously by Monnet et al. who found an unacceptable PE of 51% [22]. Second, we found that Nexfin<sup>®</sup> is most accurate in the subgroup of patients with a high

Table 2. → Comparison of several neurological, respiratory and hemodynamic variables and dose of the used drugs between the start and the end (at 8 hours) of the study period. CVP: central venous pressure, EVLWI: extravascular lung water index, GCS: Glasgow coma scale, GEDVI: global enddiastolic volume index, GEF: global ejection fraction, IAP: intra-abdominal pressure, MAP: mean arterial pressure, PPV: pulse pressure variation, SAS: sedation and agitation scale, SVRI: systemic vascular resistance index, SVV: stroke volume variation, TDCO: thermodilution cardiac output.

CO and low SVRI, however it was least accurate in patients with low CO and high SVRI. In contrast to other uncalibrated monitoring devices, NexCO keeps comparable reliability in unstable patients with severe hypotension and in patients with reduced vessel compliance due to high dose norepinephrine. In septic patients with the well known inverse TDCO/SVR hemodynamic profile, with severe hypotension or on high dose norepinephrine the calculated PE was below 30% and all criteria for interchangeability with CCO were met. Since determination of CO was only possible with invasive monitoring in the past, CO is not a target for goal directed therapy guidelines for septic patients [29]. However, since tissue perfusion and oxygen delivery are determined directly by CO and only indirectly by MAP, we strongly believe that with the development and future fine-tuning of non-invasive CO measurement devices like Nexfin®, determination of target therapy guidelines for CO should be considered [19]. Third, we found an acceptable concordance (89.3%) between the direction of changes in TDCO and NexCO during the same time interval. This was also shown in previous studies [23]. Although this analysis is based on only 3 time points and 2 values of  $\Delta$ TDCO/ $\Delta$ NexCO within an 8 hours time interval, this might be an indication that real-time measurement of the hemodynamic effects of ongoing therapeutic interventions may be reliable with the Nexfin® device. Since the absolute  $\Delta$ TDCO/ $\Delta$ NexCO amplitude did not show sufficient clinical correlation, clinical decision making should be based on hemodynamic trends rather than on absolute values of changes in measured CO. However ideally concordance should be above 90%, which was not the case when looking at changes in CCO and NexCO during the same time interval (4 values of  $\Delta$ CCO/ $\Delta$ NexCO within each 8 hours interval) with a concordance of only 81%. Fourth, analysis with polar plots showed an acceptable trending capability with 89% of the data points lying within the  $\pm 1.0$  L/min (or  $\pm 10\%$ ) limits of agreement lines. Our study is the first to use polar plot analysis in this setting. Finally, the short time to first measurement, the limited number of repositions needed to start measuring and the high scores on the nurse questionnaires illustrate that the device is very easy to use ('plug and play') in the majority of patients. However to play the devil's advocate one could also state that we were unable to use the device in 2 out of 47 patients (4.3%). Monnet et al. found worse results with the inability to obtain CO values in 15.6% of study patients [22].

	n (%)	Start	8 hours	p
<b>Neurological</b>				
Propofol (mg/kg/h)	35 (71%)	2.3±0.9	2.3±1.0	0.99
Midazolam (mg/kg/h)	31 (69%)	0.2±0.1	0.2±0.1	0.81
Remifentanyl (µg/kg/min)	39 (87%)	0.2±0.1	0.2±0.1	0.80
Cisatracurium (mg/kg/h)	9 (20%)	0.1±0.1	0.1±0.05	0.77
SAS (1—7)	45 (100%)	1.9±1.1	1.6±1.1	0.28
GCS (3—15)	45 (100%)	4.5±3.1	4.8±3.7	0.67
<b>Respiratory System</b>				
pO <sub>2</sub> /FIO <sub>2</sub>	45 (100%)	290±171	274±149	0.64
Minute ventilation (l/min)	45 (100%)	11.3±3.5	11.2±4.0	0.92
pH	45 (100%)	7.3±0.1	7.3±0.1	0.45
EVLWI (ml/kg)	45 (100%)	10.2±3.7	10.3±3.4	0.89
<b>Hemodynamics</b>				
Norepinephrine (µg/kg/min)	35 (78%)	0.2±0.2	0.2±0.2	0.90
Dobutamine (µg/kg/min)	27 (60%)	4.1±1.8	4.5±2.2	0.49
MAP (mmHg)	45 (100%)	79.6±16.9	82.9±21.1	0.41
Heart rate (BPM)	45 (100%)	90.3±25.4	88.0±23.7	0.66
CVP (mmHg)	45 (100%)	10.5±5.3	10.1±3.7	0.65
TDCO (l/min)	45 (100%)	6.2±2.1	6.9±2.2	0.14
GEF (%)	45 (100%)	22.8±8.2	24.2±8.6	0.45
GEDVI (ml/BSA)	45 (100%)	715±192	737±162	0.56
SVRI (dyne.s.cm <sup>-5</sup> /m <sup>2</sup> )	45 (100%)	1868±764	1877±793	0.56
SVV (%)	45 (100%)	15.0±8.8	12.3±7.3	0.12
PPV (%)	45 (100%)	14.6±9.7	12.1±7.9	0.17
<b>Other</b>				
IAP (mmHg)	45 (100%)	7.9±2.9	8.8±3.4	0.15
Body temperature (°C)	45 (100%)	35.4±1.7	35.5±1.7	0.76

Based on these results and review of the literature, we think the Nexfin® device can be applied in ICU or ER patients, potentially also on the regular wards and even out of hospital (if the manufacturer would

Table 3. Comparison of TDCO vs NexCO including subgroup analysis for patients with low MAP, low and high TDCO, low and high SVRI, high dose norepinephrine and hypothermia. TDCO: thermodilution cardiac output, R<sup>2</sup>: Pearson correlation coefficient, LA: limits of agreement, MAP: mean arterial pressure, SVRI: systemic vascular resistance index. \*dyne.s.cm<sup>-5</sup>/m<sup>2</sup>, NOR: norepinephrine, HYP: hypothermia.

	Unit	Overall	Low MAP	Low TDCO	High TDCO	Low SVRI	High SVRI	NOR	HYP
			≤70 mmHg	≤4 l/min	≥8 l/min	≤1700*	≥3000*	≥0,3 µg/kg/min	≤35°C
Patients	no.	45	19	9	18	29	6	11	22
Paired measurements	no.	135	27	16	32	58	14	27	63
Mean TDCO	l/min	6.6	6.2	3.8	8.2	7.2	4.5	6.7	6.5
R <sup>2</sup>	—	0.67	0.82	0.01	0.67	0.74	0.81	0.89	0.70
<i>p</i>	—	<0.001	<0.001	NS	<0.001	<0.001	<0.001	<0.001	<0.001
Bias	l/min	0.4	0.9	-0.4	0.9	0.6	-0.1	0.4	0.1
Precision	l/min	1.2	1.1	1.5	1.1	1.3	1.3	1	1.3
Lower LA	l/min	-1.9	-1.2	-3.2	-1.3	-2	-2.7	-1.5	-2.4
Upper LA	l/min	2.7	3	2.5	3	3.2	2.4	2.3	2.6
Percentage error	%	36	34	78	26	37	57	29	39

Table 4. Comparison of CCO vs NexCO including subgroup analysis for patients with low MAP, low and high TDCO, low and high SVRI, high dose norepinephrine or hypothermia. TDCO: thermodilution cardiac output, R<sup>2</sup>: Pearson correlation coefficient, LA: limits of agreement, MAP: mean arterial pressure, SVRI: systemic vascular resistance index. \*dyne.s.cm<sup>-5</sup>/m<sup>2</sup>, NOR: norepinephrine, HYP: hypothermia.

	Unit	Overall	Low MAP	Low TDCO	High TDCO	Low SVRI	High SVRI	NOR	HYP
			≤70 mmHg	≤4 l/min	≥8 l/min	≤1700*	≥3000*	≥0,3 µg/kg/min	≤35°C
Patients	no.	45	19	19	18	29	6	11	22
Paired measurements	no.	135	50	50	51	104	22	44	63
Mean TDCO	l/min	6.4	5.9	5.9	7.8	7.0	4.5	6.4	6.3
R <sup>2</sup>	—	0.71	0.87	0.87	0.73	0.81	0.78	0.94	0.67
<i>p</i>	—	<0.001	<0.001	<0.001	<0.001	<0.001	<0.001	<0.001	<0.001
Bias	l/min	0.2	0.4	0.4	0.6	0.4	0.0	0.0	0.0
Precision	l/min	1.2	0.9	0.9	1.1	1.0	1.2	0.5	1.3
Lower LA	l/min	-2.1	-1.4	-1.4	-1.7	-1.6	-2.4	-1.1	-2.7
Upper LA	l/min	2.5	2.2	2.2	2.8	2.3	2.3	1.0	2.6
Percentage error	%	36	30	30	29	29	53	16	42

provide a battery), for an initial quick hemodynamic assessment as a bridge to installation of a more advanced invasive monitoring system. Differentiation of the different types of shock on a clinical basis showed to be a major challenge and often inaccurate even in hands of experienced ICU and ER physicians [23]. Also, Nexfin® can be applied when catheter placement is problematic for instance in patients with active catheter infections after removal of the previously infected arterial line.

However, in our opinion there are 5 main reasons why Nexfin® cannot always be used as a first choice in the general ICU population with good IV and IA access requiring prolonged advanced hemodynamic monitoring. First, not only is the overall calculated PE

too high, we also think that the LA are too broad to be clinically acceptable. If a CO is measured at 8 L/min, the true value can be between 5.7 and 10.3 L/min. Of note however is that the obtained correlation coefficients and LA are comparable to previous validation studies with PiCCO against PAC [3, 15—17, 21, 28, 30] and better than results obtained with other non-calibrated, more invasive monitoring devices such as the Vigileo [26, 34], NiCO [26, 34] and PrAM [25]. Also, in a recent meta-analysis, none of the four tested methods achieved satisfactory agreement with bolus thermodilution within the expected 30% PE limits [27]. Therefore, questions are raised on the feasibility of the current validation criteria for uncalibrated CO devices. Second, in our study, Nexfin® showed to be less reliable in patients with hypothermia and to be

completely unreliable in patients with low TDCO and high SVRI (e.g. cardiogenic, obstructive and hypovolemic shock). In this subgroup NexCO showed systematic overestimation of CO. This is in line with a previous study conducted in post-cardiac surgery ICU patients showing a PE of 50%, mainly driven by inaccuracy in patients with a low CI [12]. This is in contrast however with another study in patients during CABG where half of the patients had a CI < 2.5 L/min/m<sup>2</sup> and good CO correlation coefficients and PE were still found [5]. Third, Nexfin<sup>®</sup> cannot entirely replace (less) invasive monitoring with an arterial line since arterial blood gas analyses and follow-up of lactate will always be one of the cornerstones of critical care management. Fourth, in some unstable patients and especially those with changing conditions of preload, afterload or contractility it may be advisable to calibrate the CO device in relation to the new hemodynamic situation [8]. Finally, we could not obtain any measurements in 4.3% of patients and others found that the Nexfin could not record the arterial curve due to finger hypoperfusion in 15.6% of patients [22].

Some limitations of this study need to be considered. First, this is a validation study of CO by Nexfin<sup>®</sup> against PiCCO. PAC is still considered by some clinicians as a golden standard. Although highly validated and widely used, PiCCO obtained CO shows some error against PAC. Second, the patient sample size and the size of the subgroups are probably too small to allow extensive further subgroup analysis. Third, we did not perform therapeutic intervention to assess the trending capabilities of the Nexfin. Fourth, we need to be aware that these results were obtained in an ICU patient group already receiving a lot of vasopressors and inotropic hemodynamic support thereby possibly not representing the initial hemodynamic pattern. Future studies should be performed to confirm that these results could be extrapolated to ER patients. Finally, the number of patients was determined on a random basis and no power analysis was performed.

## Conclusions

In conclusion, Nexfin<sup>®</sup> is a totally non-invasive, easy to use blood pressure and CO monitor based on finger arterial blood pressure pulse contour analysis. Nexfin<sup>®</sup> obtained CO showed a moderate to good correlation with CO measured by PiCCO although the PE was too high. The Nexfin<sup>®</sup> can be used for keeping track of changes in CO over time (for instance to assess the therapeutic effect of a given treatment), although the absolute criteria for full interchangeability were not met in this population of mixed ICU patients.

## Key messages

Nexfin<sup>®</sup> is a totally non-invasive, easy to use blood pressure and cardiac output monitor. Nexfin<sup>®</sup> shows a moderate to good cardiac output correlation with transcardiopulmonary thermodilution (TDCO) and continuous pulse contour CO (CCO) obtained by PiCCO in a mixed ICU population although the obtained percentage error was too high to allow full interchangeability. Changes in NexCO correlate well with changes in TDCO and CCO although the obtained concordance coefficient was too high to allow full trending interchangeability.

## Conflicts of interest

Manu Malbrain and Berthold Bein are members of the medical advisory board of Pulsion Medical Systems (Munich, Germany) but have no direct financial relation with the commercial identities mentioned in the paper that might lead to a conflict of interest. The other authors have no competing interests.

## Abbreviations

APACHE II: acute physiology and chronic health evaluation, CCO: Continuous cardiac output by PiCCO, CVP: central venous pressure, EVLWI: extravascular lung water index, GCS: Glasgow coma scale, GEDVI: global end diastolic volume index, GEF: global ejection fraction, IAP: intra-abdominal pressure, MAP: mean arterial pressure, NEXCO: Nexfin<sup>®</sup> obtained cardiac output, PE: percentage error, PPV: pulse pressure variation, SAPS II: simplified acute physiology score, SAS: sedation and agitation scale, SOFA: sequential organ failure assessment, SVRI: systemic vascular resistance index, SVV: stroke volume variation, TDCO: thermodilution cardiac output

**Acknowledgements** Part of this publication has been presented previously in abstract and poster form at the 31st International Symposium on Intensive Care and Emergency Medicine (ISICEM) held in Brussels in 2011 [34]. The authors are indebted to Bert Ferdinande, Nelie Desie and Annelies Gerits for their excellent input and feedback. This study won the first prize amongst the submitted abstracts on the 2nd International Fluid Academy Day, held in Antwerp, Belgium on November 17th, 2012. This study is accepted for publication in The Scientific World Journal as an open access paper under CC Licence, doi: to be defined (<http://mts.hindawi.com/author/519080.pdf>).



## References

1. Belda FJ, Aguilar G, Teboul JL, et al. PICS Investigators Group. Complications related to less-invasive haemodynamic monitoring. *Br J Anaesth* 2011;106(4):482—486
2. Biancofiore G, Critchley LA, Lee A, et al. Evaluation of an uncalibrated arterial pulse contour cardiac output monitoring system in cirrhotic patients undergoing liver surgery. *Br J Anaesth* 2009;102(1): 47—54
3. Bindels AJ, van der Hoeven JG, Graafland AD, et al. Relationships between volume and pressure measurements and stroke volume in critically ill patients. *Crit Care* 2000;4(3):193—199
4. Bland JM, Altman DG. Statistical methods for assessing agreement between two methods of clinical measurement. *Lancet* 1986;158476:307—310
5. Broch O, Renner J, Gruenewald M, et al. A comparison of the Nexfin and transcardiopulmonary thermodilution to estimate cardiac output during coronary artery surgery. *Anaesthesia* 2012; 67(4):377—383
6. Calvin JE, Driedger AA, Sibbald WJ. Does the pulmonary capillary wedge pressure predict left ventricular preload in critically ill patients? *Crit Care Med* 1981;9(6):437—443
7. Cecconi M, Rhodes A, Poloniecki J, et al. Bench-to bedside review: The importance of the precision of the reference technique in method comparison studies – with specific reference to the measurement of cardiac output. *Crit Care* 2009;13:201—206
8. Cecconi M, Malbrain ML. Cardiac output obtained by pulse pressure analysis: to calibrate or not to calibrate may not be the only question when used properly. *Intensive Care Med* 2013;39(4):787—924
9. Chen G, Meng L, Alexander B, et al. Comparison of noninvasive cardiac output measurements using the Nexfin monitoring device and the esophageal Doppler. *J Clin Anesth* 2012;24(4):275—283
10. Critchley LA, Critchley JA. A meta-analysis of studies using bias and precision statistics to compare cardiac output measurement techniques. *J Clin Monit Comput* 1999;15(2):85—91
11. Critchley LA, Lee A, Ho AM. A critical review of the ability of continuous cardiac output monitors to measure trends in cardiac output. *Anesth Analg* 2010;111(5):1180—1192
12. Fischer MO, Avram R, Cârjaliu I, et al. Non-invasive continuous arterial pressure and cardiac index monitoring with Nexfin after cardiac surgery. *Br J Anaesth* 2012;109(4):514—521
13. Genahr A, McLuckie A. Transpulmonary thermodilution in the critically ill. *Br J Int Care* 2004;1:6—10
14. Giustiniano E, Morenghez E, Ruggier N. Flotrac/Vigileo validation trials: are there reliable conclusions? *Clin Rev Opinions* 2011;3(6):62—65
15. Goepfert M, Reuter D, Akyol D, et al. Goal directed fluid management reduces vasopressor and catecholamine use in cardiac surgery patients. *Intensive Care Medicine* 2007;33:96—103
16. Holm C, Melcer B, Horbrand F, et al. Arterial thermodilution: an alternative to pulmonary artery catheter for cardiac output assessment in burn patients. *Burns* 2001;27(2):161—166
17. Kortgen A, Niederprün P, Bauer M. Implementation of an evidence-based „standard operating procedure“ and outcome in septic shock. *Crit Care Med* 2006;34(4):943—949
18. Malbrain M, De Potter P, Deeren D. Cost effectiveness of minimally invasive hemodynamic monitoring. In: Vincent J-L (Ed.) *Yearbook of Intensive Care and Emergency Medicine*, Springer Verlag, Berlin 2005; pp. 603—618
19. Malbrain MLNG, Perel A, Segal E. How to implement the recent Surviving Sepsis Campaign Guidelines at the bedside? A focus on initial fluid resuscitation. *Fluids* 2013 2(1):25—29
20. Michard F. Do we need to know cardiac preload? In: Vincent J-L (ed.) *Yearbook of Intensive Care and Emergency Medicine*, Springer Verlag, Berlin, 2004, pp. 694—701
21. Mielck F, Buhre W, Hanekop G, et al. Comparison of continuous cardiac output measurements in patients after cardiac surgery. *J Cardiothorac Vasc Anesth* 2003;17(2):211—216
22. Monnet X, Picard F, Lidzorski E, et al. The estimation of cardiac output by the Nexfin device is of poor reliability for tracking the effects of a fluid challenge. *Crit Care* 2012;16(5):R212
23. Nowak RM, Sen A, Garcia AJ, Wilkie H, Yang JJ, Nowak MR, Moyer ML. Noninvasive continuous or intermittent blood pressure and heart rate patient monitoring in the ED. *Am J Emerg Med* 2011;29(7):782—789
24. Oren-Grinberg A. The PiCCO Monitor. *Int Anesthesiol Clin* 2010;48(1):57—85
25. Paarmann H, Groesdonk H, Sedemund-Adib B, et al. Lack of agreement between pulmonary arterial thermodilution cardiac output and the pressure recording analytical method in postoperative cardiac surgery patients. *Br J Anaesth* 2011;106(4):475—481
26. Palmers PJ, Vidts W, Ameloot K, et al. Assessment of three minimally invasive continuous cardiac output measurement methods in critically ill patients and a review of the literature. *Anesthesiol Intensive Ther* 2012;44(4):188-199
27. Peyton PJ, Chong SW. Minimally invasive measurement of cardiac output during surgery and critical care. *Anesthesiology* 2010;113:220—235
28. Ritter S, Rudiger A, Maggiorini M. Transpulmonary thermodilution derived cardiac function index identifies cardiac dysfunction in acute heart failure and septic patients: an observational study. *Crit Care* 2009;13(4):R133
29. Rivers E, Nguyen B, Havstad S. Early goal-directed therapy in the treatment of severe sepsis and septic shock. *NEJM* 2001;345(19):1368—1377
30. Scheer B, Perel A, Pfeiffer UJ. Clinical review: complications and risk factors of peripheral arterial catheters used for haemodynamic monitoring in anaesthesia and intensive care medicine. *Crit Care* 2002;6(3):199—204
31. Stover JF, Stocker R, Lenherr R, et al. Noninvasive cardiac output and blood pressure monitoring cannot replace an invasive monitoring system in critically ill patients. *BMC Anesthesiology* 2009;9:6
32. Van de Vijver K, Verstraeten A, Gillebert C, et al. Validation of non-invasive hemodynamic monitoring with Nexfin in critically ill patients. *Crit Care* 2011;15(suppl. 1):P75
33. Wesseling KH, Jansen JR, Settels JJ, et al. Computation of aortic flow from pressure in humans using a nonlinear, three-element model. *J Appl Physiol* 74(5):2566—2573
34. Young BP, Low LL. Noninvasive monitoring cardiac output using partial CO(2) rebreathing. *Crit Care Clin* 2010;26(2):383—392